

(E2609*) Custom Seat Cushion Order Form

*US Only. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement.


Fields or sections with asterisk (*) are required

PO#: *	<input type="text"/>		
Client: *	<input type="text"/>		
Branch: *	<input type="text"/>	Date: *	<input type="text"/>
Contact: *	<input type="text"/>		
E-mail: *	<input type="text"/>		
Phone: *	<input type="text"/>	Fax:	<input type="text"/>
Mark for: *	<input type="text"/>	Type: *	<input type="checkbox"/> Quote <input type="checkbox"/> Order


Mobility Base Information

Manufacturer: *	<input type="text"/>	Frame Color: *	<input type="text"/>
Model Number:	<input type="text"/>	Frame Width: *	<input type="text"/>
Serial Number:	<input type="text"/>	Frame Depth: *	<input type="text"/>
		Tubing: *	<input type="checkbox"/> 3/4in (1.91cm) <input type="checkbox"/> 7/8in (2.22cm) <input type="checkbox"/> 1in (2.54cm)

CLIENT SEAT AREA MEASUREMENTS



Seat Cushion Expected Overall Dimension



F. Overall Depth: *	<input type="text"/>
G. Overall Width: *	<input type="text"/>
H. Overall Height: *	<input type="text"/>

A. Seat Width:	<input type="text"/>	C. Seat to Footplate:	<input type="text"/>
B. Knee Adduction:	<input type="text"/>	D. Total Seat Depth:	<input type="text"/>
		E. Seat to PSIS:	<input type="text"/>

1

Cushion Size

Fields or sections with asterisk (*) are required

SEAT CUSHION SIZING*

Width	10in (25.40cm)	12in (30.48cm)	14in (35.56cm)	16in (40.64cm)	18in (45.72cm)	20in (50.80cm)	22in (55.88cm)	24in (60.96cm)
Depth	10in	<input type="checkbox"/> SMD1010	<input type="checkbox"/> SMD1210					
	12in	<input type="checkbox"/> SMD1012	<input type="checkbox"/> SMD1212	<input type="checkbox"/> SMD1412				
	14in	<input type="checkbox"/> SMD1014	<input type="checkbox"/> SMD1214	<input type="checkbox"/> SMD1414	<input type="checkbox"/> SMD1614			
	16in		<input type="checkbox"/> SMD1216	<input type="checkbox"/> SMD1416	<input type="checkbox"/> SMD1616	<input type="checkbox"/> SMD1816		
	18in			<input type="checkbox"/> SMD1418	<input type="checkbox"/> SMD1618	<input type="checkbox"/> SMD1818	<input type="checkbox"/> SMD2018	
	20in				<input type="checkbox"/> SMD1620	<input type="checkbox"/> SMD1820	<input type="checkbox"/> SMD2020	<input type="checkbox"/> SMD2220
	22in					<input type="checkbox"/> SMD1822	<input type="checkbox"/> SMD2022	<input type="checkbox"/> SMD2222
	24in						<input type="checkbox"/> SMD2024	<input type="checkbox"/> SMD2224

CLIENT SPECIFIED SIZING

The seat cushion will require extended lead time if its sizing falls outside the above matrix. Set sizes by 1in (2.54cm) increments.

☐ SMDCSHSP0 Specified Dimensions Seat Base: W x D

COVER CONTACT MATERIAL*

☐ Coolcore® ☐ Air Mesh ☐ Stealth-Tek

Coolcore is a research proven performance fabric technology that adapts to its environment and regulates temperature by managing heat and moisture. The proprietary construction helps manage and distribute moisture quickly throughout the fabric to enhance the evaporation and drying time.

Air Flow Mesh, better known as **Air Mesh**, is equipped with fine pores for superior air circulation, it can assist with positioning if the client has sliding issues with the current cushion cover material. This fabric is ideal for those in hot or humid environments.

Stealth-Tek is a cast coated fabric, which are nylon/polyester fabric coated with polyurethane through the cast coating method. These type of fabrics are waterproof with stretch capability and applications in medical and healthcare. Stealth-Tek is Fluid Proof PU Film laminated to 100% polyester knit T-jersey. It is used in a variety of applications. Waterproof, Breathable and Washable.

☐ Add Polyurethane Liner (Incontinent Liner) ☐ Add Zipper Pouch

Need Help? Call +1 (512) 715-9995
Call toll free at +1 (800) 965-9229

2

Cushion Foam

FOAM OPTIONS

Lower Layer*		Middle Layer		Upper Layer	
Thickness		Thickness		Thickness	
<input type="checkbox"/> 1/2in (1.27cm)	<input type="checkbox"/> 1in (2.54 cm)	<input type="checkbox"/> 1/2in (1.27cm)	<input type="checkbox"/> 1in (2.54 cm)	<input type="checkbox"/> 1/2in (1.27cm)	<input type="checkbox"/> 1in (2.54 cm)
Foam Type		Foam Type		Foam Type	
<input type="checkbox"/> SunMate® FRG	<input type="checkbox"/> Open Cell	<input type="checkbox"/> SunMate® FRG	<input type="checkbox"/> Open Cell	<input type="checkbox"/> SunMate® FRG	<input type="checkbox"/> Open Cell
<input type="checkbox"/> Medium Soft	<input type="checkbox"/> Medium	<input type="checkbox"/> X-Soft	<input type="checkbox"/> Medium	<input type="checkbox"/> XX-Soft	<input type="checkbox"/> Medium
<input type="checkbox"/> Medium		<input type="checkbox"/> Soft		<input type="checkbox"/> X-Soft	
		<input type="checkbox"/> Medium Soft		<input type="checkbox"/> Soft	

Note: The cushion requires a 5/8in (1.59cm) closed cell foam base which adds to the preceding selections. The cover included with the cushion supports up to 4in (10.16cm) of foam thickness, any thickness beyond that will require extended lead time.




3 Basic Cushion Shape

Fields or sections with asterisk (*) are required

3.1 Select Base Format*

Use increments of 1in (2.54cm)



☐ **Anti-Thrust (AT)**

Anti-Thrust Depth

Anti-Thrust Height

☐ 1/2in (1.27cm)


☐ 1in (2.54cm)

☐ 1 1/2in (3.81cm)

☐ 2in (5.08cm)

☐ 2 1/2in (6.35cm)

OR



☐ **Omit Anti-Thrust (Planar)**


A. Cushion Back Height:

B. Cushion Front Height:

Apply only to Anti-Thrust. Specify total cushion front height expected.


OR

3.2 Select Contour Option



☐ **Anti-Thrust with Adductor**

Expected adduction total height



Adductor Specifications

C Adductor Height (from Surface)


☐ 1/2in (1.27cm)

☐ 1in (2.54cm)


☐ 1 1/2in (3.81cm)

☐ 2in (5.08cm)

☐ 2 1/2in (6.35cm)




☐ **Planar with Adductor**



☐ **Anti-Thrust with Abductor**

Expected Abduction total height



Abductor Specifications


D Abductor Depth
Set at 1in (2.54cm) increments

E Abductor Height (from Surface)
Increments of 1/2in (1.27cm)


☐ 1in (2.54cm)

☐ 1 1/2in (3.81cm)

☐ 2in (5.08cm)




☐ **Planar with Abductor**



☐ **Full Contour with Anti-Thrust**

Contoured Seat Cushion

Full contoured seat cushion requires A, C, D and E field values. If you selected Anti-Thrust it will require B field value as well.



☐ **Planar Full Contour**

Need Help? Call +1 (512) 715-9995
Call toll free at +1 (800) 965-9229

** For reference only (not required)

3.1 Special Cut Outs

LEG LENGTH DISCREPANCY

☐ **Short Right** (select depth below)

☐ **Short Left** (select depth below)

☐ 1/2in (1.27cm)
 ☐ 1in (2.54cm)
 ☐ 1 1/2in (3.81cm)
 ☐ 2in (5.08cm)

GROWTH ADAPTATIONS

☐ **I Style**

☐ **Growth Notch** (Select Depth)

☐ 1in (2.54cm)
 ☐ 2in (5.08cm)
 ☐ 3in (7.62cm)

RELIEFS

☐ **Hamstring Relief**

☐ **Ischial Relief** (select gel type below)


☐ Viscous Gel
 ☐ Ultra-Lite Gel

Good for End Users who have poor trunk control and need better immersion to stabilize the pelvis.
 Good for end users who need more movement and less shear. Pediatric users for example.


4

Mounting

- ☐ No Mounting
- ☐ Seat Pan




Seat Pan Hardware ☐ No Hardware



- ☐ Trouble Free
- ☐ Knob
- ☐ J & L

Mounting Format: ☐ 4 pt ☐ 2 pt

Optional:



☐ Seat-to-Back

Order the bracket when 2 point mounting format is selected, and the current back system does not include one. Also, consider if the seating system expects no flexing by the posterior area.

Offset: ☐ 1 1/2in (3.81cm) ☐ 3in (7.62cm) ☐ Adjustable




Also available
ADI™
Carbon Fiber
Seat Pans
for weights
under 250 lb



2 Point Mounting comes with seat tabs for posterior support.

☐ Rigidizer ACM ABS Contour Mold



- ☐ **Small**
16in Width (40.64cm)
- ☐ **Medium**
18in Width (45.72cm)
- ☐ **Large**
20in Width (50.80cm)

Not compatible with the I-Style Cut Out and might not be compatible with some cushion widths or depths, see note below.

Can be cut on designated ridges to the following sizes:

Small: 16in x 16in, 16in x 18in
Medium: 18in x 16in, 18in x 18in, 18in x 20in
Large: 20in x 18in, 20in x 20in

Need Help? Call +1 (512) 715-9995
Call toll free at +1 (800) 965-9229

6 Part Number Reference

COVERS / FOAMS / GELS

Covers

SMDCOV	10	10	C
Cushion Width	10 to 24		
Cushion Depth	10 to 24		
Material:			
C: Coolcore®			
M: Airmesh			
T: Stealth-Tek			

Polyurethane Liners

SMDPUL	10	10
Cushion Width	10 to 24	
Cushion Depth	10 to 24	

Foam Options

SMDFO	#	-	#	-	#	-	#
Width Option:							
N: 10in to 12in							
S: 13in to 16in							
M: 17in to 20in							
L: 21in to 24in							
Layer Position:							
U: Upper (Top)							
M: Middle							
L: Lower (Bottom)							
Foam Type:							
OC3: Open Cell Medium							
SMFRG1: SunMate® FRG Extra Extra Soft							
SMFRG2: SunMate FRG Extra Soft							
SMFRG3: SunMate FRG Soft							
SMFRG4: SunMate FRG Medium Soft							
SMFRG5: SunMate FRG Medium							
Thickness:							
0.5: 1/2in (1.27cm)							
1: 1in (2.54cm)							

IT Gel

SMDINS	GUL	S
Gel Type:		
GUL: Gel Ultra Lite		
VIS: Viscous Gel		
Size (based on width)		
S: Small (10in to 16in)		
L: Large (16in to 24in)		

STRUCTURE

Anti-Thrust

SMDANT	N	06	-	0.5
Size Width:				
N: 10in to 12in				
S: 13in to 16in				
M: 17in to 20in				
L: 21in to 24in				
Size Depth:				
06: up to 6in (15.24cm)				
08: up to 8in (20.32cm)				
10: up to 10in (25.40cm)				
12: up to 12in (30.48cm)				
Height:				
0.5: 1/2in (1.27cm)				
1.0: 1in (2.54cm)				
1.5: 1 1/2in (3.81cm)				
2.0: 2in (5.08cm)				
2.5: 2 1/2in (6.35cm)				

Abductor

SMDABD	06	-	1
Size Depth:			
06: up to 6in (15.24cm)			
08: up to 8in (20.32cm)			
10: up to 10in (25.40cm)			
12: up to 12in (30.48cm)			
Height:			
1.0: 1in (2.54cm)			
1.5: 1 1/2in (3.81cm)			
2.0: 2in (5.08cm)			

Hamstring Relief

SMDHSR	1
Options:	
1: 1in (2.54cm) Depth 45° Slope	
2: 2in (5.08cm) Depth 45° Slope	

Adductor

SMDADD	N	0.5
Length:		
N: 10in to 12in		
S: 13in to 16in		
M: 17in to 20in		
L: 21in to 24in		
Height:		
0.5: 1/2in (1.27cm)		
1.0: 1in (2.54cm)		
1.5: 1 1/2in (3.81cm)		
2.0: 2in (5.08cm)		
2.5: 2 1/2in (6.35cm)		

Leg Length Discrepancy

SMDLLD	0.5	-	L
Depth:			
0.5: 1/2in (1.27cm)			
1.0: 1in (2.54cm)			
1.5: 1 1/2in (3.81cm)			
2.0: 2in (5.08cm)			
2.5: 2 1/2in (6.35cm)			
Side:			
L: Left			
R: Right			

Growth Notch

SMDGN	1
Options:	
1: 1in (2.54cm) Depth	
2: 2in (5.08cm) Depth	

MOUNTING

Rigidizer

SMDRGD	SML
Options:	
SML: 16in (40.64cm) Wide	
MED: 18in (45.72cm) Wide	
LRG: 20in (50.80cm) Wide	






Seat Pans

SMDSP	10	10
Cushion Width	10 to 24	
Cushion Depth	10 to 24	

Mounting Hardware

SMDHW	JL	2	P	-	034
TYPE:					
JL: J & L					
KN: Knob					
TF: Trouble Free					
Set Point:					
2: 2 Point					
4: 4 Point					
Tubing:					
034: for 3/4in (1.91cm) Tubing					
078: for 7/8in (2.22cm) Tubing					
100: for 1in (2.54cm) Tubing					

CUSTOM CUSHION SKETCH



This image shows a full page of blank, lined paper. It features approximately 28 horizontal blue or grey lines spaced evenly apart, typical of notebook paper. The lines extend across the entire width of the page, leaving small margins at the top and bottom. There are no vertical lines, text, or other markings present.



+1 (512) 715-9995
+1 (800) 965-9229 Toll Free
+1 (512) 715-9954 Fax
+1 (800) 806-1225 Toll Free Fax

104 John Kelly Drive,
Burnet, TX 78611
information@sw.stealthproducts.com
www.stealthproducts.com
facebook.com/StealthProducts
twitter.com/stealthproducts