





STANDARD CONFIGURATIONS



The UniLink lateral system is a versatile and highly adjustable support solution, adaptable to the varied and changing requirements of users. UniLink thoracic adductors are available in two standard lengths, long or short, and can be configured as fixed or swing-away, providing increased accessibility for patients and attendants. Extra links are available to extend the reach of the laterals if needed, and a custom pad size option helps ensure the perfect fit. Mounting options are available for a variety of powered and manual mobility systems.




Description	Part Number	A	B	C	D
Short Fixed Lateral (L/R)	ULS23200-(L/R)	1.0" (25mm)	2.5" (63mm)	3.25" (80 mm)	1.0" (25.4mm)
Long Fixed Lateral (L/R)	ULS23210-(L/R)			5" (126mm)	
Short Swing-Away Lateral (L/R)	ULS25000-(L/R)			4.5" (115mm)	
Swing-Away Lateral (L/R)	ULS25100-(L/R)			6.25" (159mm)	
In-Line Swing-Away Lateral (L/R)	ULS05100-(L/R)			6" (152mm)	
Short Fixed Lateral (L/R) FSC Option*	ULS2320F-(L/R)		3.0" (78mm)	3.25" (80mm)	
Long Fixed Lateral (L/R) FSC Option*	ULS2321F-(L/R)			5" (126 mm)	
Short Swing-Away Lateral (L/R) FSC Option*	ULS2500F-(L/R)			4.5" (115mm)	
Swing-Away Lateral (L/R) FSC Option*	ULS2510F-(L/R)			6.25" (159mm)	
In-Line Swing Away Lateral (L/R) FSC Option*	ULS0510F-(L/R)			6" (152mm)	

* FSC Option requires additional part U1PFS1100 listed in "Extra Parts" on the reverse of this page.



Stealth Tech Dashboard

UniLink™ Thoracic Lateral Hardware



*UIMBKIP3G
Permobil Corpus 3G
Bracket Kit
(two included in kit)*

*UIMBK1048
Permobil 3G
Bracket Only*


*UIMBK1090
90° Bracket*

MOUNTING OPTIONS


Quantum TB3 Kit	<i>UIMBKITB3</i>
Quantum TB4 Kit	<i>UIMBKITB4</i>
TB3/TB4 Bracket only	<i>UIMBK1125</i>
Permobil 3G Kit	<i>UIMBKIP3G</i>
Permobil 3G Bracket Only	<i>UIMBK1048</i>
90° Bracket (Planar/TFB I-Back, Diego, Flex Back)	<i>UIMBK1090</i>
105° Bracket (TFB R-Back)	<i>UIMBK1105</i>
135° Bracket (ADI, Q-Back, TFB Contour)	<i>UIMBK1135</i>
3/4" Cane Clamp	<i>UIMCNI340</i>
7/8" Cane Clamp	<i>UIMCNI780</i>
1" Cane Clamp	<i>UIMCNI100</i>
3/4" Crescent Clamp*	<i>UIMCRI340</i>
7/8" Crescent Clamp*	<i>UIMCRI780</i>
1" Crescent Clamp*	<i>UIMCRI100</i>
All Season Bracket - 70°	<i>UIMASI070</i>
All Season Bracket - 75°	<i>UIMASI075</i>
All Season Bracket - 80°	<i>UIMASI080</i>
All Season Bracket - 90°	<i>UIMASI090</i>



*UIMCRI100
Crescent Clamp
* currently unavailable*




*UIMASI090
All Season
Bracket*



*UIMCNI100
Cane
Clamp*

EXAMPLE CONFIGURATIONS



*All Season
mounting.*

TECHNICAL DATA

All UniLink fasteners are metric,
and most can be installed or adjusted
with a 4mm Hex Key.




ULP535050320

LATERAL PAD OPTIONS


3x5	<i>ULP535050320</i>
3x5 Gel	<i>ULP535057120</i>
4x6	<i>ULP546050320</i>
4x6 Gel	<i>ULP546057120</i>
5x7	<i>ULP557050320</i>
5x7 Gel	<i>ULP557057120</i>
Custom	<i>ULPCUS10110</i>

EXTRA PARTS


Additional 1.5" Link	<i>UILNK150</i>
Additional 2" Link	<i>UILNK200</i>
Additional 3" Link	<i>UILNK300</i>
Additional 2" In-Line Swing-Away Link	<i>UILIS200</i>
Lateral Full Surface Contact (required for hardware packages ULS2320F, ULS2321F, ULS2500F, ULS2510F & ULS0510F)	<i>UIPFSI100</i>



UILNK150



UILIS200



UIPFSI100