





STANDARD CONFIGURATIONS




The UniLink lateral system is a versatile and highly adjustable support solution, adaptable to the varied and changing requirements of users. UniLink Thigh Guides are available in long or short models, have several mounting options, and can be configured as fixed or swing-away. Extra links are available to extend the reach of the laterals if needed, and a custom pad size option helps ensure the perfect fit. Mounting options are available for a variety of powered and manual mobility systems.



Description	Part Number	A	B	C	D	E
Short Stem Thigh Guide Track Package	UTTP5-Q	1" (2.5mm)	1.77" (45 mm)	4.53" (115 mm)	5" (127 mm)	N/A
Long Stem Thigh Guide Track Package	UTTP7-Q		3" (76 mm)	6.57" (167 mm)		
Adjustable Thigh Guide (L/R)	UTT2311902-(L/R)		N/A	N/A	N/A	5.5" (140 mm)
Short Multi-Angle Thigh Guide (L/R)	UTT7A11602-(L/R)		1.77" (45 mm)	4.53" (115 mm)		
Short Stave Thigh Guide (L/R)	UTT7A1260000				5.38" (137 mm)	N/A
Long Multi-Angle Thigh Guide (L/R)	UTT6A11602-(L/R)		3" (76 mm)	6.57" (167 mm)		
Long Stave Thigh Guide (L/R)	UTT6A1260000					
Short Flip-Away Thigh Guide (L/R)*	UTS25100-(L/R)		N/A	N/A	N/A	7" (178 mm)
Long Flip-Away Thigh Guide (L/R)*	UTS25300-(L/R)					8.5" (216 mm)

* FSC available only on flip-away models.





TRACK PAD OPTIONS

4x4	UTP544056100
4x6	UTP546056100
4x8	UTP548056100
4x10	UTP540056100
4x12	UTP54D056100
Custom*	UTP5CUS06100

* Please call for quote on custom orders above 4x12



NEOPRENE PAD OPTIONS

3x6	UTP536050300
4x4	UTP544050300
4x6	UTP546050300
4x8	UTP548050300
4x10	UTP540050300
4x12	UTP54D050300
Custom*	UTP5CUS50300

* Please call for quote on custom orders above 4x12

EXAMPLE CONFIGURATIONS




TECHNICAL DATA


All UniLink fasteners are metric, and most can be installed or adjusted with a 4mm Hex Key.

MOUNTING OPTIONS

Pelvic/Thigh Bracket - 90°	UIMPBI090
Pelvic/Thigh Bracket - Cane Mount 3/4"	UIMPCI340
Pelvic/Thigh Bracket - Cane Mount 7/8"	UIMPCI780
Pelvic/Thigh Bracket - Cane Mount 1"	UIMPCI100
TB4 Receiver	UT800000-Q
Pelvic/Thigh Bracket - All Season	UTMASI090




UTPFSI100
Full Surface
Contact Option




UIMPCI100
Cane
Clamp

EXTRA PARTS

Additional 1.5" Link	UILNK150
Additional 2" Link	UILNK200
Additional 3" Link	UILNK300
Additional 2" In-Line Swing-Away Link	UILIS200
Lateral Full Surface Contact (only available on UTS25100 & UTS25300)	UTPFSI100



UILIS200
Inline Latch



UILNK150
1.5" Link