







## STANDARD CONFIGURATIONS



Quantum TB4 armrest  
with  
**UETB4021FR-Q**  
(TB4 Bracket & Lateral)  
and  
**UEP530050320-R**  
(3x10 Elbow Pad)



*The UniLink lateral system is a versatile and highly adjustable support solution, adaptable to the varied and changing requirements of users. UniLink's Elbow Block is intended to provide extra support of the elbow and forearm by keeping the arm steady and the hand in range of the control device during tilt actions or maneuvers over rough terrain. The blocks come in fixed or swing-away configurations with a range of pad sizes to meet the needs of users and technicians. Mounting options are available for a variety of powered mobility systems.*



Description	Part Number	A
Fixed Elbow Block*	UE000021F0	4" (100 mm)
Swing-Away Elbow Block* (Left, Right)	UE000213S(L,R)	3.5" (89 mm)
TB4 Standard Fixed Elbow Block (Left, Right)	UETB4021F(L,R)-Q	4" (100 mm)
TB4 Standard Swing-Away Elbow Block (Left, Right)	UETB4213S(L,R)-Q	3.5" (89 mm)
Swing-Away Elbow Block with 3x3 Pad* (Left, Right)	UEPKG213S(L,R)-33	3.5" (89 mm)
Swing-Away Elbow Block with 4x4 Pad* (Left, Right)	UEPKG213S(L,R)-44	3.5" (89 mm)
Swing-Away Elbow Block with 3x10 Pad* (Left, Right)	UEPKG213S(L,R)-30	3.5" (89 mm)

\* Hardware kit does not include a mount. Please see Mounting Options on the reverse of this sheet.



1" mount with Permobil nuts  
UEPRMTRNFO

### MOUNTING OPTIONS

1" Mount with Permobil Nuts (Left/Right)	UEPRMTRNF(L,R)
1" Mount with Single Track Nuts (Left/Right)	UESTDTRNF(L,R)



Swing-Away Elbow Block  
UE0000213S-L

### ELBOW PADS

3x3 Elbow Pad, Neoprene	UEP533050320
3x3 Elbow Pad, Gel	UEP533057120
4x4 Elbow Pad, Neoprene	UEP544050320
4x4 Elbow Pad, Gel	UEP544057120
3x10 Elbow Pad, Neoprene, Left	UEP53005032L
3x10 Elbow Pad, Gel, Left	UEP53005712L
3x10 Elbow Pad, Neoprene, Right	UEP53005032R
3x10 Elbow Pad, Gel, Right	UEP53005712R

### EXAMPLE CONFIGURATIONS

1" Mount with  
Permobil Nuts  
(UEPRMTRNFL),  
Fixed Elbow Block  
Hardware (UE000021FL)  
and 3x3 Neoprene Pad  
(UEP533050320)




1" Mount with  
Permobil Nuts  
(UEPRMTRNFL),  
Fixed Elbow Block  
Hardware (UE000021FL)  
and 3x3 Neoprene Pad  
(UEP544050320)



### TECHNICAL DATA

All UniLink fasteners are metric,  
and most can be installed or adjusted  
with a 4mm Hex Key.



### ACCESSORIES

*Stealth's Armrest Side Support is an excellent adjunct to the UniLink Elbow Block Hardware. In conjunction with a 3x3 or 4x4 pad, the Armrest Side Support provides six inches of side support for the arm on the armrest.*



### ARMREST SIDE SUPPORT

Armrest Side Support 6 x 2 CoolCore	ASF6X2C-S-Q
-------------------------------------	-------------

