



CUSTOM SEATING SOLUTIONS





SOLUTIONS FOR LIFE

Stealth Product's seating offerings have evolved in quality, innovation, and turn time since its inception over 12 years ago, now making Stealth the premier provider of custom seating solutions in the industry.



STEALTH CUSTOM SEAT [E2609]

A custom seat solution that is fully configurable, from structural modifications, to our proprietary FRG foams and Coolcore fabric, Stealth custom seats can be configured to meet the specific level of support and positioning required by the end-user. All of Stealth's custom seat and back solutions enable you to incorporate Stealth secondary supports for additional postural needs.




STEALTH TRUE FORMING BACK [E2617]



A custom back solution designed for ease of use, configurability, and client function.

Starting with client specific dimensions, Stealth custom TFB Backs are configured to meet the exact postural needs of the individual. In addition, Stealth incorporates unique core materials into our custom seating solutions to enhance comfort/sitting tolerance, as well as elements designed to manage surface climate.

Ease of use is accomplished by minimizing additional hardware required to incorporate secondary positioning components, such as lateral thoracic supports. The unique shell design enables providers to direct mount positioning components to the shell of the back via universal adjustment slots and pre-drilled attaching holes. This results in an easy to install, low maintenance solution, that will function as well as it looks.



TARTA ORIGINAL [E2617]

In addition to our lineup of True Forming Back custom solutions, Stealth also offers a unique customizable solution with the Tarta Original back. This modular design allows for virtually limitless configurability:

- Height adjustable, Flexible spine
- Posterior and Lateral flex in the back
- Range of Malleable Ribs to accommodate size and contour
- Wide selection of pads

This is a great alternative to some of the custom molding being utilized to address fixed asymmetrical positioning. Tarta is lighter weight and adaptable to change, as well as configurable on-site for immediate delivery.



104 John Kelly Drive
Burnet, Texas 78611

(800) 965-9229

stealthproducts.com

