



Stealth Glacial SP Cushion Order Form

104 John Kelly Drive, Burnet, TX 78611 | Phone: 800-965-9229 | Fax: 800-806-1225 | Email: info@stealthproducts.com

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Marked for: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are on-chair MSRP. Prices, specifications, part numbers, and availability are subject to change without notice.

GLACIAL SP CUSHION

This cushion offers skin protection and a 4-way stretch, breathable cover. Weight capacities are 175 lbs for 10"-14" width, 350 lbs for 15"-21" width, and 400 lbs for 22" width and greater. The XYYY portion of all part numbers on this page directly relate to the dimensions of the cushion, where XX represents the width and YY represents the depth. For example, a 16" wide by 18" deep cushion would be ST-GLSP1618. Available sizes are shown on the chart to the right.

1. Cushion Selection

- Glacial SP Cushion, 10"-20" Wide, E2603 \$347
- Glacial SP Cushion, 22"-24" Wide, E2604 \$462

Cushion Part: **ST-GLSP**XXYY

Cushion Width: _____

Cushion Depth: _____

| | | Width | | | | | | | | | | | |
|-------|-----|-------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| | | 10" | 12" | 14" | 15" | 16" | 17" | 18" | 19" | 20" | 21" | 22" | 24" |
| Depth | 10" | | | | - | - | - | - | - | - | - | - | - |
| | 12" | | | | - | - | - | - | - | - | - | - | - |
| | 14" | | | | - | - | - | - | - | - | - | - | - |
| | 15" | - | - | - | | - | - | - | - | - | - | - | - |
| | 16" | - | - | | | | - | | | | | | - |
| | 17" | - | - | | | | | | | | | | - |
| | 18" | - | - | | | | | | | | | | |
| | 19" | - | - | | | | | | | | | | |
| | 20" | - | - | | | | | | | | | | |
| | 21" | - | - | | | | | | | | | | |
| | 22" | - | - | | | | | | | | | | |
| 23" | - | - | | | | | | | | | | | |
| 24" | - | - | | | | | | | | | | | |

2. Add Fluid Proof Liner

The Glacial SP Cushion is available with an optional Fluid Proof Urethane or Fluid Proof Polyurethane Liner.

2.1a. Install Fluid Proof Urethane Liner

- Credit for Urethane Liner (for 21"W and less) -\$53
Part: PTOINDV3927
- Credit for Urethane Liner (22"W and greater) -\$63
Part: PTOINDV3928

2.1b. Fluid Proof Urethane Liner

| Cushion Size | Part | MSRP | Cushion Size | Part | MSRP | Cushion Size | Part | MSRP |
|----------------|----------------------|------|----------------|----------------------|------|----------------|----------------------|------|
| 10"W x 10"D... | URL1709070-1010..... | \$53 | 17"W x 17"D... | URL1709070-1717..... | \$53 | 20"W x 22"D... | URL1709070-2022..... | \$63 |
| 10"W x 12"D... | URL1709070-1012..... | \$53 | 17"W x 18"D... | URL1709070-1718..... | \$53 | 21"W x 19"D... | URL1709070-2119..... | \$63 |
| 10"W x 14"D... | URL1709070-1014..... | \$53 | 17"W x 19"D... | URL1709070-1719..... | \$53 | 21"W x 20"D... | URL1709070-2120..... | \$63 |
| 12"W x 10"D... | URL1709070-1210..... | \$53 | 17"W x 20"D... | URL1709070-1720..... | \$53 | 21"W x 21"D... | URL1709070-2121..... | \$63 |
| 12"W x 12"D... | URL1709070-1212..... | \$53 | 18"W x 16"D... | URL1709070-1816..... | \$53 | 22"W x 16"D... | URL1709070-2216..... | \$63 |
| 12"W x 14"D... | URL1709070-1214..... | \$53 | 18"W x 17"D... | URL1709070-1817..... | \$53 | 22"W x 18"D... | URL1709070-2218..... | \$63 |
| 14"W x 10"D... | URL1709070-1410..... | \$53 | 18"W x 18"D... | URL1709070-1818..... | \$53 | 22"W x 20"D... | URL1709070-2220..... | \$63 |
| 14"W x 12"D... | URL1709070-1412..... | \$53 | 18"W x 20"D... | URL1709070-1820..... | \$53 | 22"W x 22"D... | URL1709070-2222..... | \$63 |
| 14"W x 14"D... | URL1709070-1414..... | \$53 | 18"W x 22"D... | URL1709070-1822..... | \$53 | 22"W x 24"D... | URL1709070-2224..... | \$63 |
| 14"W x 16"D... | URL1709070-1416..... | \$53 | 19"W x 17"D... | URL1709070-1917..... | \$53 | 24"W x 18"D... | URL1709070-2418..... | \$63 |
| 14"W x 18"D... | URL1709070-1418..... | \$53 | 19"W x 18"D... | URL1709070-1918..... | \$53 | 24"W x 19"D... | URL1709070-2419..... | \$63 |
| 15"W x 15"D... | URL1709070-1515..... | \$53 | 19"W x 19"D... | URL1709070-1919..... | \$53 | 24"W x 20"D... | URL1709070-2420..... | \$63 |
| 15"W x 16"D... | URL1709070-1516..... | \$53 | 19"W x 20"D... | URL1709070-1920..... | \$53 | 24"W x 21"D... | URL1709070-2421..... | \$63 |
| 15"W x 17"D... | URL1709070-1517..... | \$53 | 20"W x 16"D... | URL1709070-2016..... | \$53 | 24"W x 22"D... | URL1709070-2422..... | \$63 |
| 16"W x 16"D... | URL1709070-1616..... | \$53 | 20"W x 18"D... | URL1709070-2018..... | \$53 | 24"W x 23"D... | URL1709070-2423..... | \$63 |
| 16"W x 18"D... | URL1709070-1618..... | \$53 | 20"W x 20"D... | URL1709070-2020..... | \$53 | 24"W x 24"D... | URL1709070-2424..... | \$63 |
| 16"W x 20"D... | URL1709070-1620..... | \$53 | | | | | | |

Fluid Proof Liner section is continued on the following page.

2. Add Fluid Proof Liner continued

2.2a. Install Fluid Proof Polyurethane Liner

Install Polyurethane Liner \$0
 Part: PTOINDV3918

2.2b. Fluid Proof Polyurethane Liner

| Cushion Size | Part | MSRP | Cushion Size | Part | MSRP | Cushion Size | Part | MSRP |
|----------------|----------------------|------|----------------|----------------------|------|----------------|----------------------|------|
| 10"W x 10"D... | PUL1709030-1010..... | \$58 | 17"W x 17"D... | PUL1709030-1717..... | \$58 | 20"W x 22"D... | PUL1709030-2022..... | \$58 |
| 10"W x 12"D... | PUL1709030-1012..... | \$58 | 17"W x 18"D... | PUL1709030-1718..... | \$58 | 21"W x 19"D... | PUL1709030-2119..... | \$58 |
| 10"W x 14"D... | PUL1709030-1014..... | \$58 | 17"W x 19"D... | PUL1709030-1719..... | \$58 | 21"W x 20"D... | PUL1709030-2120..... | \$58 |
| 12"W x 10"D... | PUL1709030-1210..... | \$58 | 17"W x 20"D... | PUL1709030-1720..... | \$58 | 21"W x 21"D... | PUL1709030-2121..... | \$58 |
| 12"W x 12"D... | PUL1709030-1212..... | \$58 | 18"W x 16"D... | PUL1709030-1816..... | \$58 | 22"W x 16"D... | PUL1709030-2216..... | \$79 |
| 12"W x 14"D... | PUL1709030-1214..... | \$58 | 18"W x 17"D... | PUL1709030-1817..... | \$58 | 22"W x 18"D... | PUL1709030-2218..... | \$79 |
| 14"W x 10"D... | PUL1709030-1410..... | \$58 | 18"W x 18"D... | PUL1709030-1818..... | \$58 | 22"W x 20"D... | PUL1709030-2220..... | \$79 |
| 14"W x 12"D... | PUL1709030-1412..... | \$58 | 18"W x 20"D... | PUL1709030-1820..... | \$58 | 22"W x 22"D... | PUL1709030-2222..... | \$79 |
| 14"W x 14"D... | PUL1709030-1414..... | \$58 | 18"W x 22"D... | PUL1709030-1822..... | \$58 | 22"W x 24"D... | PUL1709030-2224..... | \$79 |
| 14"W x 16"D... | PUL1709030-1416..... | \$58 | 19"W x 17"D... | PUL1709030-1917..... | \$58 | 24"W x 18"D... | PUL1709030-2418..... | \$79 |
| 14"W x 18"D... | PUL1709030-1418..... | \$58 | 19"W x 18"D... | PUL1709030-1918..... | \$58 | 24"W x 19"D... | PUL1709030-2419..... | \$79 |
| 15"W x 15"D... | PUL1709030-1515..... | \$58 | 19"W x 19"D... | PUL1709030-1919..... | \$58 | 24"W x 20"D... | PUL1709030-2420..... | \$79 |
| 15"W x 16"D... | PUL1709030-1516..... | \$58 | 19"W x 20"D... | PUL1709030-1920..... | \$58 | 24"W x 21"D... | PUL1709030-2421..... | \$79 |
| 15"W x 17"D... | PUL1709030-1517..... | \$58 | 20"W x 16"D... | PUL1709030-2016..... | \$58 | 24"W x 22"D... | PUL1709030-2422..... | \$79 |
| 16"W x 16"D... | PUL1709030-1616..... | \$58 | 20"W x 18"D... | PUL1709030-2018..... | \$58 | 24"W x 23"D... | PUL1709030-2423..... | \$79 |
| 16"W x 18"D... | PUL1709030-1618..... | \$58 | 20"W x 20"D... | PUL1709030-2020..... | \$58 | 24"W x 24"D... | PUL1709030-2424..... | \$79 |
| 16"W x 20"D... | PUL1709030-1620..... | \$58 | | | | | | |

3. Alternate Cushion Cover

The Glacial SP cushion is standard with a Coolcore cover which can be substituted with a Mesh or Stealth Tek Cover. A credit will be applied.

3.1a. Install Mesh Cover

Credit for Mesh as Installed Cover-\$105
 Part: PTOINDV3887

3.1b. Mesh Cover Selection

| Cushion Size | Part | MSRP | Cushion Size | Part | MSRP | Cushion Size | Part | MSRP |
|-----------------|--------------------|-------|-----------------|--------------------|-------|-----------------|--------------------|-------|
| 10"W x 10"D.... | GLSP-MBLK1010..... | \$105 | 17"W x 17"D.... | GLSP-MBLK1717..... | \$105 | 20"W x 22"D.... | GLSP-MBLK2022..... | \$105 |
| 10"W x 12"D.... | GLSP-MBLK1012..... | \$105 | 17"W x 18"D.... | GLSP-MBLK1718..... | \$105 | 21"W x 19"D.... | GLSP-MBLK2119..... | \$105 |
| 10"W x 14"D.... | GLSP-MBLK1014..... | \$105 | 17"W x 19"D.... | GLSP-MBLK1719..... | \$105 | 21"W x 20"D.... | GLSP-MBLK2120..... | \$105 |
| 12"W x 10"D.... | GLSP-MBLK1210..... | \$105 | 17"W x 20"D.... | GLSP-MBLK1720..... | \$105 | 21"W x 21"D.... | GLSP-MBLK2121..... | \$105 |
| 12"W x 12"D.... | GLSP-MBLK1212..... | \$105 | 18"W x 16"D.... | GLSP-MBLK1816..... | \$105 | 22"W x 16"D.... | GLSP-MBLK2216..... | \$105 |
| 12"W x 14"D.... | GLSP-MBLK1214..... | \$105 | 18"W x 17"D.... | GLSP-MBLK1817..... | \$105 | 22"W x 18"D.... | GLSP-MBLK2218..... | \$105 |
| 14"W x 10"D.... | GLSP-MBLK1410..... | \$105 | 18"W x 18"D.... | GLSP-MBLK1818..... | \$105 | 22"W x 20"D.... | GLSP-MBLK2220..... | \$105 |
| 14"W x 12"D.... | GLSP-MBLK1412..... | \$105 | 18"W x 20"D.... | GLSP-MBLK1820..... | \$105 | 22"W x 22"D.... | GLSP-MBLK2222..... | \$105 |
| 14"W x 14"D.... | GLSP-MBLK1414..... | \$105 | 18"W x 22"D.... | GLSP-MBLK1822..... | \$105 | 22"W x 24"D.... | GLSP-MBLK2224..... | \$105 |
| 14"W x 16"D.... | GLSP-MBLK1416..... | \$105 | 19"W x 17"D.... | GLSP-MBLK1917..... | \$105 | 24"W x 18"D.... | GLSP-MBLK2418..... | \$105 |
| 14"W x 18"D.... | GLSP-MBLK1418..... | \$105 | 19"W x 18"D.... | GLSP-MBLK1918..... | \$105 | 24"W x 19"D.... | GLSP-MBLK2419..... | \$105 |
| 15"W x 15"D.... | GLSP-MBLK1515..... | \$105 | 19"W x 19"D.... | GLSP-MBLK1919..... | \$105 | 24"W x 20"D.... | GLSP-MBLK2420..... | \$105 |
| 15"W x 16"D.... | GLSP-MBLK1516..... | \$105 | 19"W x 20"D.... | GLSP-MBLK1920..... | \$105 | 24"W x 21"D.... | GLSP-MBLK2421..... | \$105 |
| 15"W x 17"D.... | GLSP-MBLK1517..... | \$105 | 20"W x 16"D.... | GLSP-MBLK2016..... | \$105 | 24"W x 22"D.... | GLSP-MBLK2422..... | \$105 |
| 16"W x 16"D.... | GLSP-MBLK1616..... | \$105 | 20"W x 18"D.... | GLSP-MBLK2018..... | \$105 | 24"W x 23"D.... | GLSP-MBLK2423..... | \$105 |
| 16"W x 18"D.... | GLSP-MBLK1618..... | \$105 | 20"W x 20"D.... | GLSP-MBLK2020..... | \$105 | 24"W x 24"D.... | GLSP-MBLK2424..... | \$105 |
| 16"W x 20"D.... | GLSP-MBLK1620..... | \$105 | | | | | | |

3.2a. Install Stealth Tek Cover

Credit for Stealth Tek as Alternate Cover-\$105
 Part: PTOINDV3886

Alternate Tek 2 Cover section is continued on the following page.

3. Alternate Cushion Cover continued

3.2b. Stealth Tek Cover Selection

| Cushion Size | Part | MSRP | Cushion Size | Part | MSRP | Cushion Size | Part | MSRP |
|-----------------|---------------------|-------|-----------------|---------------------|-------|-----------------|---------------------|-------|
| 10"W x 10"D.... | GLSP-TBLK1010 | \$105 | 17"W x 17"D.... | GLSP-TBLK1717 | \$105 | 21"W x 19"D.... | GLSP-TBLK2119 | \$105 |
| 10"W x 12"D.... | GLSP-TBLK1012 | \$105 | 17"W x 18"D.... | GLSP-TBLK1718 | \$105 | 21"W x 20"D.... | GLSP-TBLK2120 | \$105 |
| 10"W x 14"D.... | GLSP-TBLK1014 | \$105 | 17"W x 19"D.... | GLSP-TBLK1719 | \$105 | 21"W x 21"D.... | GLSP-TBLK2121 | \$105 |
| 12"W x 10"D.... | GLSP-TBLK1210 | \$105 | 17"W x 20"D.... | GLSP-TBLK1720 | \$105 | 22"W x 16"D.... | GLSP-TBLK2216 | \$105 |
| 12"W x 12"D.... | GLSP-TBLK1212 | \$105 | 18"W x 16"D.... | GLSP-TBLK1816 | \$105 | 22"W x 18"D.... | GLSP-TBLK2218 | \$105 |
| 12"W x 14"D.... | GLSP-TBLK1214 | \$105 | 18"W x 17"D.... | GLSP-TBLK1817 | \$105 | 22"W x 20"D.... | GLSP-TBLK2220..... | \$105 |
| 14"W x 10"D.... | GLSP-TBLK1410 | \$105 | 18"W x 18"D.... | GLSP-TBLK1818 | \$105 | 22"W x 22"D.... | GLSP-TBLK2222..... | \$105 |
| 14"W x 12"D.... | GLSP-TBLK1412 | \$105 | 18"W x 20"D.... | GLSP-TBLK1820 | \$105 | 22"W x 24"D.... | GLSP-TBLK2224..... | \$105 |
| 14"W x 14"D.... | GLSP-TBLK1414 | \$105 | 18"W x 22"D.... | GLSP-TBLK1822 | \$105 | 24"W x 18"D.... | GLSP-TBLK2418 | \$105 |
| 14"W x 16"D.... | GLSP-TBLK1416 | \$105 | 19"W x 17"D.... | GLSP-TBLK1917 | \$105 | 24"W x 19"D.... | GLSP-TBLK2419 | \$105 |
| 14"W x 18"D.... | GLSP-TBLK1418 | \$105 | 19"W x 18"D.... | GLSP-TBLK1918 | \$105 | 24"W x 20"D.... | GLSP-TBLK2420..... | \$105 |
| 15"W x 15"D.... | GLSP-TBLK1515 | \$105 | 19"W x 19"D.... | GLSP-TBLK1919 | \$105 | 24"W x 21"D.... | GLSP-TBLK2421 | \$105 |
| 15"W x 16"D.... | GLSP-TBLK1516 | \$105 | 19"W x 20"D.... | GLSP-TBLK1920 | \$105 | 24"W x 22"D.... | GLSP-TBLK2422..... | \$105 |
| 15"W x 17"D.... | GLSP-TBLK1517 | \$105 | 20"W x 16"D.... | GLSP-TBLK2016..... | \$105 | 24"W x 23"D.... | GLSP-TBLK2423 | \$105 |
| 16"W x 16"D.... | GLSP-TBLK1616 | \$105 | 20"W x 18"D.... | GLSP-TBLK2018..... | \$105 | 24"W x 24"D.... | GLSP-TBLK2424..... | \$105 |
| 16"W x 18"D.... | GLSP-TBLK1618 | \$105 | 20"W x 20"D.... | GLSP-TBLK2020..... | \$105 | | | |
| 16"W x 20"D.... | GLSP-TBLK1620 | \$105 | 20"W x 22"D.... | GLSP-TBLK2022..... | \$105 | | | |

4. Add Solid Base

The Glacial SP Cushion is available with an optional ABS Flat Rigidizer or ACM Rigidizer (For use with manual wheelchair only).

4.1. ABS Flat Rigidizer

| Cushion Size | Part | MSRP | Cushion Size | Part | MSRP | Cushion Size | Part | MSRP |
|----------------|-----------------------|------|----------------|-----------------------|------|----------------|--------------------|------|
| 10"W x 10"D... | DWR1010E188A..... | \$48 | 17"W x 17"D... | DWR1010E234H | \$48 | 20"W x 16"D... | DWR1010E193G..... | \$58 |
| 10"W x 12"D... | DWR1010E188C..... | \$48 | 17"W x 18"D... | DWR1010E234H | \$48 | 20"W x 18"D... | DWR1010E193I..... | \$58 |
| 10"W x 14"D... | DWR1010E188E..... | \$48 | 17"W x 19"D... | DWR1010E234J..... | \$48 | 20"W x 20"D... | DWR1010E193K..... | \$58 |
| 12"W x 10"D... | DWR1010E189A..... | \$48 | 17"W x 20"D... | DWR1010E234H | \$48 | 20"W x 22"D... | DWR1010E193M..... | \$58 |
| 12"W x 12"D... | DWR1010E189C..... | \$48 | 18"W x 16"D... | DWR1010E192G..... | \$48 | 22"W x 16"D... | DWR1010E194G..... | \$58 |
| 12"W x 14"D... | DWR1010E189E..... | \$48 | 18"W x 17"D... | DWR1010E192G..... | \$48 | 22"W x 18"D... | DWR1010E194I..... | \$58 |
| 14"W x 10"D... | MEC1804152-1410 | \$48 | 18"W x 18"D... | DWR1010E192I..... | \$48 | 22"W x 20"D... | DWR1010E194K..... | \$58 |
| 14"W x 12"D... | DWR1010E190C..... | \$48 | 18"W x 20"D... | DWR1010E192K..... | \$48 | 22"W x 22"D... | DWR1010E194M | \$58 |
| 14"W x 14"D... | DWR1010E190E..... | \$48 | 18"W x 22"D... | DWR1010E192M | \$48 | 22"W x 24"D... | DWR1010E193O..... | \$58 |
| 14"W x 16"D... | DWR1010E190G..... | \$48 | 19"W x 17"D... | MEC1804152-1917 | \$48 | 24"W x 18"D... | DWR1010E236I | \$58 |
| 14"W x 18"D... | DWR1010E190I..... | \$48 | 19"W x 18"D... | MEC1804152-1917 | \$48 | 24"W x 19"D... | DWR1010E236I | \$58 |
| 15"W x 15"D... | MEC1804152-1515 | \$48 | 19"W x 19"D... | MEC1804152-1919 | \$48 | 24"W x 20"D... | DWR1010E236K..... | \$58 |
| 15"W x 16"D... | MEC1804152-1515 | \$48 | 19"W x 20"D... | MEC1804152-1919 | \$48 | 24"W x 21"D... | DWR1010E236K..... | \$58 |
| 15"W x 17"D... | DWR1010E232H | \$48 | 21"W x 19"D... | MEC1804152-2119 | \$48 | 24"W x 22"D... | DWR1010E236M | \$58 |
| 16"W x 16"D... | DWR1010E191G..... | \$48 | 21"W x 20"D... | MEC1804152-2119 | \$48 | 24"W x 23"D... | DWR1010E236M | \$58 |
| 16"W x 18"D... | DWR1010E191I..... | \$48 | 21"W x 21"D... | MEC1804152-2121 | \$48 | 24"W x 24"D... | DWR1010E236O | \$58 |
| 16"W x 20"D... | DWR1010E191K..... | \$48 | | | | | | |

4.2. ACM Rigidizer

For use with manual wheelchair only.

| Cushion Size | Part | MSRP | Cushion Size | Part | MSRP | Cushion Size | Part | MSRP |
|-----------------|----------------------|-------|-----------------|----------------------|-------|-----------------|----------------------|-------|
| 16"W x 16"D ... | CUS142467-1616 | \$100 | 17"W x 20"D... | CUS142466 | \$100 | 19"W x 18"D ... | CUS142465-2018 | \$100 |
| 16"W x 18"D ... | CUS142467..... | \$100 | 18"W x 16"D ... | CUS142466-1816 | \$100 | 19"W x 19"D ... | CUS142465-2018 | \$100 |
| 16"W x 20"D... | CUS142467..... | \$100 | 18"W x 17"D... | CUS142466-1816 | \$100 | 19"W x 20"D... | CUS142465 | \$100 |
| 17"W x 17"D... | CUS142466-1816 | \$100 | 18"W x 18"D... | CUS142466-1818 | \$100 | 20"W x 18"D... | CUS142465-2018 | \$100 |
| 17"W x 18"D... | CUS142466-1818 | \$100 | 18"W x 20"D... | CUS142466 | \$100 | 20"W x 20"D... | CUS142465 | \$100 |
| 17"W x 19"D... | CUS142466-1818 | \$100 | 18"W x 22"D... | CUS142466 | \$100 | 20"W x 22"D... | CUS142465 | \$100 |