



Measurement and Fitment Reference



**CONTOUR
BACK**



**ACTIVE
DEEP
EXTRA-DEEP**

** Extra Deep Option shown **

SHELL SPECIFICATIONS BASED ON WIDTH

(All measurements in inches)

| Shell Dimension Width | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 |
|----------------------------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|
| Active Shell Width (A) | 7.4 | 8.4 | 9.3 | 10.3 | 11.1 | 12.0 | 13.0 | 14.0 | 15.0 | 16.0 | 17.0 | 18.0 | 19.0 | 20.0 | 21.0 |
| 10-14W Wing Width (B) | 7.5 | 8.5 | 9.4 | 10.4 | 11.2 | 12.2 | 13.2 | 14.2 | 15.2 | 16.2 | 17.2 | 18.2 | 19.2 | 20.2 | 21.2 |
| 15-22W Wing Width (B) | 7.5 | 8.5 | 9.4 | 10.4 | 12.6 | 13.6 | 14.6 | 15.6 | 16.6 | 17.6 | 18.6 | 19.6 | 20.6 | 21.6 | 22.6 |
| Active Depth (C) | 1.8 | 2.3 | | 2.6 | | | | | | | 2.7 | | | | |
| Deep Depth (C) | 3.8 | 4.3 | | 4.6 | | | | | | | 4.7 | | | | |
| Ex. Deep Depth (C) | 5.8 | 6.3 | | 6.6 | | | | | | | 6.7 | | | | |
| Active Lat. Wd. (D) | 6.7 | 7.7 | 8.6 | 9.6 | 10.4 | 11.3 | 12.3 | 13.3 | 14.3 | 15.3 | 16.3 | 17.3 | 18.3 | 19.3 | 20.3 |
| 10-14W Dp/XDp Lat. Wd. (D) | 5.5 | 6.5 | 7.4 | 8.4 | 9.2 | 10.2 | 11.2 | 12.2 | 13.2 | 14.2 | 15.2 | 16.2 | 17.2 | 18.2 | 19.2 |
| 15-22W Dp/XDp Lat. Wd. (D) | 5.5 | 6.5 | 7.4 | 8.4 | 10.6 | 11.6 | 12.6 | 13.6 | 14.6 | 15.6 | 16.6 | 17.6 | 18.6 | 19.6 | 20.6 |
| Hardware Min. Width | 10.3 | 11.3 | 12.2 | 13.0 | 14.0 | 15.0 | 16.0 | 17.0 | 18.0 | 19.0 | 20.0 | 21.0 | 22.0 | 23.0 | 24.0 |
| Hardware Max. Width | 13.1 | 14.2 | 15.0 | 16.0 | 16.8 | 17.8 | 18.8 | 19.8 | 20.8 | 21.8 | 22.8 | 23.8 | 24.8 | 25.8 | 26.8 |

- Height dimensions on Contour shells are true-to-size, i.e. TFBMD1818CCA is 18 in. tall. Overall foam height will be 1 inch larger than the shell.
- Figures shown for (D) in the table above are calculated with 1.5 in. foam. Other foam thicknesses may alter the lateral width specification.