

+1 (512) 715-9995
+1 (800) 965-9229 Toll Free
+1 (512) 715-9954 Fax
+1 (800) 806-1225 Toll Free Fax

104 John Kelly Drive,
Burnet, TX 78611
information@sw.stealthproducts.com
www.stealthproducts.com
facebook.com/StealthProducts
twitter.com/stealthproducts

This order uses Adobe® Acrobat® interactive technology and we recommend you to use latest version of Acrobat Reader after downloading and saving the form to your desktop local drive.

US Only. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement.

[E2617] TFB107 - True Forming Linear Back

Client Information

PO#:¹

Client:*

Branch:*

Date:

Contact:*

E-mail:²

Phone:²

Fax:

Mark for:*

Order Type:* ☐ Quote ☐ Order

Mobility Base Information

Manufacturer:

Model:

Serial:

Frame Color:

Frame Width:

Frame Depth:

Cane Size:* ☐ 3/4in (1.91cm) ☐ 7/8in (2.22cm) ☐ 1in (2.54cm)

Mount Set Up:* ☐ Between Canes ☐ Full Width

*Required Fields
1. Required if this is an order
2. One of the options is required

1 Back Support Size*

	Back Width							
	10	12	14	16	18	20	22	24
22	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
18	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
14	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
12	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
10	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Price increase will show from 21in (55.88cm) width and/or height and above. Anything above 24in (60.96cm) width or 22in (55.88cm) height or under 10in (25.40cm) either width or height call a representative.


☐ Other: Width: Height:

Price increase will show from 21in (55.88cm) width and/or height and above. Anything above 24in (60.96cm) width or 22in (55.88cm) height or under 10in (25.40cm) either width or height call a representative.

1.1 Back Format*

Back Shell & Cushion
Reference Number
TFBMD

☐ R-Back



2 Cover Choices*

☐ Coolcore® ☐ StealthTek

☐ Airmesh

All covers will have Honeycomb¹ for visible non-contact fabric and Loop² material on the back of cushion.

3 Cover Options

☐ Privacy Flap TFBMD051 \$97.85
Connected with zipper attachment to the back bottom cover for privacy.

☐ Embroidery \$0.00

Thread Color: ☐ Blue ☐ Pink ☐ Red ☐ Green ☐ Purple ☐ Silver

Font Style: ☐ Block ☐ Script

Embroidery Name:

Embroidery TFBMDOM

4 Custom Foam Options*

Predefined Foam Configuration

Contact your Customer Representative when ordering a Predefined Foam Selection.

Base Layer

☐ HR Medium \$0 ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

SunMate® FRG Foams \$180.25

☐ Extra Extra Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Extra Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Medium Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Medium ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

Middle Layer

optional

☐ HR Medium \$0 ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

SunMate® FRG Foams \$180.25


☐ Extra Extra Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Extra Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Medium ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Medium Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)



Contact Layer

optional

☐ HR Medium \$0 ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

SunMate® FRG Foams \$180.25

☐ Extra Extra Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Extra Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Medium Soft ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

☐ Medium ☐ 1/2in (1.27cm) ☐ 1in (2.54cm)

Foam Structure Breakdown

Contact Foam TFBF - C -

Middle Foam TFBF - M -

Base Foam TFBF - B -

5 Foam-in-Place

Total Back Thickness	Frame Width	MSRP
<input type="checkbox"/> 3in (7.62cm)	14in to 16in (35.56cm to 40.64cm)	\$303.85
<input type="checkbox"/> 3in (7.62cm)	Over 20in (50.80cm)	\$607.70

† Refers to expected total back thickness.

Includes:

1. Back cushion cover with optional 3in (7.62cm) to 4in (10.16cm) thickness with elastic pleated surface for foam growth
2. 1in (2.54cm) Medium FRG SunMate® base foam layer
3. 1/4in (1.27cm) Extra Soft SunMate contact foam layer
4. 1 plastic bag for use during molding process

F.I.P.S is setup as middle layer foam

6 Back Ergonomics

- ☐ Lumbar Support Insert TFBMDOL10 \$139.05
- ☐ Curved by Foam TFBMDOF12 \$105.00

Need Help? Call +1 (512) 715-9995
Call toll free at +1 (800) 965-9229

6.2 Bi-Angular Back Options

One Piece Bi-Angular Back
TFBMDOB1 \$303.85
Split shell, solid cushion and cover

☐ 3in (7.62cm)
☐ 4in (10.16cm)
☐ 5in (12.70cm)
☐ 6in (15.24cm)
☐ 7in (17.78cm)

Two Piece Bi-Angular Back
TFBMDOB2 \$365.65
Two cushion pieces covered independent from each other

☐ 3in (7.62cm)
☐ 4in (10.16cm)
☐ 5in (12.70cm)
☐ 6in (15.24cm)
☐ 7in (17.78cm)

7 Lateral Support

Pad Options \$105.45 ea. Custom \$168.71

E0956

Depth

Height

Lateral Pad(s) Part Number
TFB- -L
TFB- -R

For custom sizes above 6in (12.70cm) in height or 8in (20.32cm) in depth please call for a quote.

Pad Size (Height x Depth)

<input type="checkbox"/> 3in x 5in (7.62cm x 12.70cm)	<input type="checkbox"/> Left <input type="checkbox"/> Right
<input type="checkbox"/> 4in x 6in (12.70cm x 15.24cm)	<input type="checkbox"/> Left <input type="checkbox"/> Right
<input type="checkbox"/> 5in x 7in (12.70cm x 17.78cm)	<input type="checkbox"/> Left <input type="checkbox"/> Right
<input type="checkbox"/> Custom: Height x Depth	<input type="checkbox"/> Left <input type="checkbox"/> Right

Pad Style:

☐ Flat ☐ Curved

Laterals Cover Options*

☐ Coolcore® ☐ StealthTek

Laterals come with Fixed Mounting Bracket to attach it to back.

6.1 R-Back Wings

Lower Wings K0108

☐ Left **TFBMDLW 0**

☐ Extra Small ☐ Large
☐ Small ☐ Extra Large
☐ Medium

☐ Right **TFBMDLW 1**

☐ Extra Small ☐ Large
☐ Small ☐ Extra Large
☐ Medium

[1] Extra Small	10in (25.40cm) height and over	\$180.25
[2] Small	14in (35.56cm) height and over	\$190.55
[3] Medium	16in (40.64cm) height and over	\$200.85

Bi-Angular cut out availability is dependent of back height see table below for reference

B. Bi-Angular Height:	3in (7.62cm)	4in (10.16cm)	5in (12.70cm)	6in (15.24cm)	7in (17.78cm)
A. Back Height					
10 in (25.40cm)	✓	✓			
11in (27.24cm)	✓	✓			
12in (30.48cm)	✓	✓	✓		
13in (33.02cm)	✓	✓	✓		
14in (35.56cm)	✓	✓	✓	✓	
15in (38.10cm)	✓	✓	✓	✓	
16in (40.64cm)	✓	✓	✓	✓	✓
17in (43.18cm)	✓	✓	✓	✓	✓
18in (45.72cm)	✓	✓	✓	✓	✓
19in (48.26cm)	✓	✓	✓	✓	✓
20in (50.80cm) and over	✓	✓	✓	✓	✓

✓ Available

7.1 Laterals Mounting Solutions

Mounting Options E1028

☐ TWB Link Style \$316.21 ea.

☐ Left ☐ Right

☐ Thin \$316.21 ea.

☐ Left ☐ Right ☐ w/ Standard Link ☐ w/ Standard Link

☐ Modular \$350.20 ea.

☐ Left ☐ Right

Lateral Hardware Breakdown
TFB- -L
TFB- -R

Lateral Solutions Notes:

R-Back Wings include pad. Pad cover materials same as back. ☐ Roll Up

Upper Wings K0108

☐ Left **TFBMDLW 2**

☐ Extra Small ☐ Large
☐ Small ☐ Extra Large
☐ Medium

☐ Right **TFBMDLW 3**

☐ Extra Small ☐ Large
☐ Small ☐ Extra Large
☐ Medium

[4] Large	18in (45.72cm) height and over	\$211.15
[5] Extra Large	20in (50.80cm) height and over	\$221.45

Bi-Angular cut out availability is dependent of back height see table below for reference

B. Bi-Angular Height:	3in (7.62cm)	4in (10.16cm)	5in (12.70cm)	6in (15.24cm)	7in (17.78cm)
A. Back Height					
10 in (25.40cm)	✓	✓			
11in (27.24cm)	✓	✓			
12in (30.48cm)	✓	✓	✓		
13in (33.02cm)	✓	✓	✓		
14in (35.56cm)	✓	✓	✓	✓	
15in (38.10cm)	✓	✓	✓	✓	
16in (40.64cm)	✓	✓	✓	✓	✓
17in (43.18cm)	✓	✓	✓	✓	✓
18in (45.72cm)	✓	✓	✓	✓	✓
19in (48.26cm)	✓	✓	✓	✓	✓
20in (50.80cm) and over	✓	✓	✓	✓	✓

✓ Available

☐ Roll Up Roll up lateral pads and hardware

☐ UniLink Swing-Away Hdw \$316.21 ea.

☐ Left ☐ Right

☐ UniLink In-Line Swing-Away Hdw \$307 ea.

☐ Left ☐ Right

☐ UniLink Swing-Away Hdw FSC \$307 ea.

Requires FSC Mount

☐ Left ☐ Right

☐ UniLink Swing-Away Hdw FSC Short \$307 ea.

Requires FSC Mount

☐ Left ☐ Right

☐ TFB Style, Fixed Hdw \$119.48 ea.

☐ Left ☐ Right

☐ TFB Style, UniLink Swing-Away Hdw Short \$307 ea.

☐ Left ☐ Right

☐ Lateral Full Surface (FSC) Mount

Requires one for each FSC Hdw selected above.

HCPCS:K0108 \$108.88 ea.

8 Back Mounting Options*

J&L Hardware

☐ TFBHWRJL2 2 Points \$53.05

☐ TFBHWRJL4 4 Points \$53.05

Knob Hardware (All Purpose)

☐ TFBHWRAK2 2 Points \$297.05

☐ TFBHWRAK4 4 Points \$328.88

Mounting Notes:

9 Head Solution

Head Pad CP CCFPM

Comfort Plus

☐ 6in (15.24cm) \$170.25

6 1/4in (15.87cm) Width x 4 1/4in (10.76cm) Height

☐ 8in (20.32cm) \$170.25

9 9/16in (24.25cm) Width x 5 3/4in (14.60cm) Height

☐ 10in (25.40cm) \$170.25

11 7/32in (28.51cm) Width x 5 3/4in (14.60cm) Height

☐ 14in (35.56cm) \$170.25

14 5/32in (35.94cm) Width x 5 3/4in (14.60cm) Height

☐ 19in (48.26cm) \$259.96

19 1/2in (49.53cm) Width x 5 1/8in (13.01cm) Height

☐ Custom: Width X Height \$350.09

For custom headrest pads over 20in (50.80cm) Width and 5.75in (14.61cm) Height call representative for a quote.

Head Solutions Notes:

Trouble Free (TF)

☐ TFBHWRTF2 2 Points \$53.05

☐ TFBHWRTF4 4 Points \$53.05

ADI Quick Release

☐ TFBHWRQR2F 2 Points \$53.05

Configuration might change depending on selections. Only available as 2 points Full Width.

9.1 Head Solution Mounting

Basic Mounting

☐ Basic \$179.10

☐ Multi-Axis Basic \$215.27

Height Adjustment

☐ 14in (35.56cm) Standard Stem

☐ 10in (25.40cm) Short Stem

Light Multi-Axis Mounting

☐ Multi-Axis \$414.76

Multi-Axis

PWC Mounting

☐ PWC Mnt; Quantum TB3/TB4 \$509.85

SSAPT B4

☐ PWC Mnt; Quantum TB3/TB4, Pediatric (14W and less) \$509.85

SSAPT B4

☐ PWC Mnt, General (Permobil/Rovi) \$509.85

SSAPGEN

☐ PWC Mnt, General (Permobil/Rovi) Pediatric (14W and less) \$509.85

SSAPGEN-PED

Compatible with Quantum Rehab® TRU-Balance® 3 or 4 and Edge® Stretto™

Power Wheelchair Direct Mounting Interface Laser Cutout for Permobil, Quantum and Sunrise PWC \$311.06

Precision laser cutouts in the back shell allow mounting directly to the power wheelchair's frame reducing seat depth loss. Available for Invacare, Permobil, Quantum and Sunrise power wheelchairs.

☐ TFBMNT-PERM Permobil

☐ TFBMNT-QNTM Quantum

☐ TFBMNT-SUNR Sunrise

Not compatible with I-Backs use SSAPT B4 option above.

☐ Roll Up Roll up hardware

Hardware Part Number TFB-HS

UniLink Style Mounting

☐ UniLink™ \$316.21

1/2in squared bracket

☐ UniLink™ Flip Down \$428.48

Height Adjustment

☐ 12in (30.48cm) Standard Stem

☐ 7in (17.78cm) Short Stem

TWB Link Style Mounting

☐ TWB Link Style with reversed depth offset \$302.82

15mm squared bracket

☐ TWB Link Style \$315.18

☐ TWB Link Style with Flip Down \$428.48

Height Adjustment

☐ 14in (35.56cm) Standard Stem

☐ 10in (25.40cm) Short Stem

Height Adjustment

☐ 10in (25.54cm) Standard Stem

☐ 6in (15.24cm) Short Stem

Height Adjustment

☐ 10in (25.54cm) Standard Stem

☐ 6in (15.24cm) Short Stem