


# Measurement and Fitment Reference




\*\* Extra Deep Option shown \*\*

**CONTOUR  
BACK**



**ACTIVE  
DEEP  
EXTRA-DEEP**



## SHELL SPECIFICATIONS BASED ON WIDTH

(All measurements in inches)

Stated Width	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
Actual Shell Width (A)	7.4	8.4	9.3	10.2	11.2	12.2	13.2	14.2	15.2	16.2	17.2	18.2	19.2	20.2	21.2	
<b>Active</b> Depth (C)	1.8	2.3		2.6	2.7											
<b>Deep</b> Depth (C)	3.8	4.3		4.6	4.7											
<b>Ex. Deep</b> Depth (C)	5.8	6.3		6.6	6.7											
<b>Active</b> Lateral Width (D)	6.7	7.7	8.6	9.6	10.4	11.3	12.3	13.3	14.3	15.3	16.3	17.3	18.3	19.3	20.3	
<b>10-14W Dp/XDp</b> Lateral Width (D)	5.5	6.5	7.4	8.4	9.2	10.2	11.2	12.2	13.2	14.2	15.2	16.2	17.2	18.2	19.2	
<b>15-22W Dp/XDp</b> Lateral Width (D)	5.5	6.5	7.4	8.4	10.6	11.6	12.6	13.6	14.6	15.6	16.6	17.6	18.6	19.6	20.6	
<b>Fixed HW</b>	Min. Width	10.2	11.2	12.2	13.0	14.0	15.0	16.0	17.0	18.0	19.0	20.0	21.0	22.0	23.0	24.0
	Max. Width	13.0	14.0	15.0	15.8	16.8	17.8	18.8	19.8	20.8	21.8	22.8	23.8	24.8	25.8	26.8
<b>Removable HW</b>	Min. Width	8.7	9.9	10.9	11.9	12.5	13.4	14.4	15.4	16.4	17.4	18.4	19.4	20.4	21.4	22.4
	Max. Width	10.7	11.9	12.9	13.9	14.5	15.4	16.4	17.4	18.4	19.4	20.4	21.4	22.4	23.4	24.4

- Height dimensions on Contour shells are true-to-size, i.e. TFBMD1818CCA is 18 in. tall. Overall foam height will be 1 inch larger than the shell.
- Figures shown for (D) in the table above are calculated with 1.5 in. foam. Other foam thicknesses may alter the lateral width specification.