



# Measurement and Fitment Reference



\*\* Active Option shown \*\*



**FLEX  
BACK**



**ACTIVE  
DEEP**

## SHELL SPECIFICATIONS BASED ON HEIGHT

(All measurements in inches)

Shell Dimension Height	10	11	12	13	14	15	16	17	18	19	20	21	22
Center Cutout Height (C)	5.5	6.5	7.5	7.5	8.5	9.0	10.0	11.0	11.0	10.5	11.5	12.5	13.5

## SHELL SPECIFICATIONS BASED ON WIDTH

BASED ON HEIGHT (All measurements in inches)

Stated Width	12	13	14	15	16	17	18	19	20	21	22	23	24	
Top Width (A)	4.0	5.0	6.0	7.0	8.0	9.0	10.0	11.0	12.0	13.0	14.0	15.0	16.0	
Max Width (B)	8.5	9.5	10.5	11.5	12.5	13.5	14.5	15.5	16.5	17.5	18.5	19.5	20.5	
Center Cutout Width (D)	2.0	3.0	4.0	5.0	5.0	6.0	7.0	8.0	9.0	10.0	11.0	12.0	13.0	
<b>Active</b> Depth (E)	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	
<b>Deep</b> Depth (E)	6.1	6.1	6.1	6.1	6.1	6.1	6.1	6.1	6.1	6.1	6.1	6.1	6.1	
<b>Active</b> Lateral Width (F)	5.6	6.6	7.6	8.6	9.6	10.6	11.6	12.6	13.6	14.6	15.6	16.6	17.6	
<b>Deep</b> Lateral Width (F)	7.8	8.8	9.8	10.8	11.8	12.8	13.8	14.8	15.8	16.8	17.8	18.8	19.8	
<b>Fixed HW</b>	Hardware Min. Width	11.5	12.5	13.5	14.5	15.5	16.5	17.5	18.5	19.5	20.5	21.5	22.5	23.5
	Hardware Max. Width	14.3	15.3	16.3	17.3	18.3	19.3	20.3	21.3	22.3	23.3	24.3	25.3	26.3
<b>Removable HW</b>	Hardware Min. Width	10.2	11.2	12.2	13.2	14.2	15.2	16.2	17.2	18.2	19.2	20.2	21.2	22.2
	Hardware Max. Width	12.2	13.2	14.2	15.2	16.2	17.2	18.2	19.2	20.2	21.2	22.2	23.2	24.2

- Height dimensions on Flex shells are true-to-size, i.e. TFBMD1818FMA is 18 in. tall. Overall foam height will be 1 inch larger than the shell.
- Figures shown for (F) in the table above are calculated with 1.5 in. foam. Other foam thicknesses may alter the lateral width specification.